



## MOUSE ANTI-KERATIN PAN (4, 5, 6, 8, 10, 13, 18) MONOCLONAL ANTIBODY

<b>CATALOG NUMBER:</b>	MAB1636
<b>LOT NUMBER:</b>	
<b>QUANTITY:</b>	100 $\mu$ L
<b>SPECIFICITY:</b>	Human keratins 4, 5, 6, 8, 10, 13, 18. Reacts with intermediate filaments of many species up to <i>Xenopus laevis</i> .
<b>ISOTYPE:</b>	IgG <sub>1</sub>
<b>APPLICATIONS:</b>	Immunocytochemistry on cultured cells (methanol/acetone fixed). Immunohistochemistry on frozen tissue. Suitable for use on paraffin embedded tissue (protease digestion may be necessary). Western blot. Optimal working dilutions must be determined by end user.
<b>FORMAT:</b>	Ascites fluid.
<b>STORAGE/HANDLING:</b>	Store at -20°C in undiluted aliquots for up to 12 months. Avoid repeated freeze/thaw cycles.
<b>REFERENCES:</b>	1. <i>Tumor Marker Oncol.</i> (1990) <b>5</b> :219.

**Important Note:** *During shipment, small volumes of product will occasionally become entrapped in the seal of the product vial. For products with volumes of 200  $\mu$ L or less, we recommend gently tapping the vial on a hard surface or briefly centrifuging the vial in a tabletop centrifuge to dislodge any liquid in the container's cap.*

FOR RESEARCH USE ONLY; NOT FOR USE IN DIAGNOSTIC  
PROCEDURES. NOT FOR HUMAN OR ANIMAL CONSUMPTION

Unless otherwise stated in our catalog or other company documentation accompanying the product(s), our products are intended for research use only and are not to be used for any other purpose, which includes but is not limited to, unauthorized commercial uses, in vitro diagnostic uses, ex vivo or in vivo therapeutic uses or any type of consumption or application to humans or animals.

©2002 - 2010: Millipore Corporation. All rights reserved. No part of these works may be reproduced in any form without permission in writing.

28820 Single Oak Drive • Temecula, CA 92590  
Technical Support: T: 1-800-MILLIPORE (1-800-645-5476) • F: 1-800-437-7502  
[www.millipore.com](http://www.millipore.com)