

Product Information

Anti-Human Serum

produced in rabbit, whole antiserum

Catalog Number **H8765**

Product Description

Anti-Human Serum was produced in rabbit using pooled normal human serum as the immunogen. The antiserum has been treated to remove lipoproteins.

Reagent

Supplied as a liquid containing 15 mM sodium azide as preservative.

Precautions and Disclaimer

This product is for R&D use only, not for drug, household, or other uses. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

Storage/Stability

For continuous use, store at 2-8 °C for up to one month. For extended storage, the solution may be frozen in working aliquots. Repeated freezing and thawing, or storage in "frost-free" freezers, is not recommended. If slight turbidity occurs upon prolonged storage, clarify the solution by centrifugation before use.

Product Profile

The protein concentration is determined by Biuret.

Strong reactivity with normal human serum has been determined by immunoelectrophoresis (IEP). This antiserum has not been assayed for interspecies cross reactivity.

DS,PHC 06/15-1