



## BRUCELLA SELECTIVE SUPPLEMENT

Product Number **B3176**

### Product Description

Brucella Selective Supplement is an antibiotic supplement recommended for the selective isolation of the *Brucella* species from milk.

### Components

#### Item

Polymyxin B	2,500 IU
Bacitracin	12,500 IU
Nystatin	50,000 IU
Cycloheximide	50.0 mg
Nalidixic Acid	2.5 mg
Vancomycin	10.0 mg

**Precautions and Disclaimer** For laboratory use only.

Not for drug, household or other uses.

## Product Information

### Preparation Instructions

Rehydrate the contents of 1 vial aseptically with 5 mls of 50% (M3641) methanol. Shake to form a uniform suspension. Each vial is sufficient for 500 mls of medium.

### Storage Temperature

Store at 2-8°C.

Sigma brand products are sold through Sigma-Aldrich, Inc.

Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.