

Product Information

Anti-SFRS6

Developed in Rabbit, Affinity isolated antibody

Product Number **S 2820**

Product Description

Anti-SFRS6 is developed in rabbit using as immunogen a synthetic peptide, RSKDEYEKRSR corresponding to human SFRS6 (amino acids 264-275) conjugated to BSA. The antibody is affinity-purified using the immunizing peptide immobilized on resin.

Anti-SFRS6 specifically recognizes SFRS6 by immunoblotting (39.6 kDa) and immunohistochemistry. The antibody detects human and rodent SFRS6. Other species reactivity has not been confirmed.

SFRS6 (Splicing factor arginine/serine-rich 6), aliases pre-mRNA splicing factor SRP55, B52 and MGC5045, is a member of the SR protein family of splicing factors. SFRS6 is a regulator of alternative RNA splicing.

Reagent

Anti-SFRS6 is provided as affinity isolated antibody in a 50% ammonium sulfate suspension in phosphate buffered saline, containing no additional preservatives.

Preparation Instructions

Method 1 for immunostaining and immunoblotting (Western blot)

1. Carefully resuspend antibody pellet to uniformity.
2. Remove a fixed amount of suspension and dissolve 1:10 in PBS or TBS to yield a 100 µg/ml solution.

Method 2 for immunoprecipitation, supershift, immunostaining and immunoblotting (Western blot)

1. Pellet antibodies at 10,000 – 15,000 x g for 10 minutes at 2 to 8 °C using a microcentrifuge.

2. Carefully remove as much supernatant as possible. It is not necessary to remove all the ammonium sulfate solution; a small residual amount will not effect the antibody preparation. Dissolve the pellet (antibody) in small volume (100 µL) of PBS (or TBS) at final concentration of 1 mg/ml (100 µg/100 µl). Do not allow the pellet to dry out. This can cause loss of activity. Gently allow pellet to dissolve at least 1 hour before use. Do not vortex. Mix by finger-tapping or gentle stirring.

Notes:

- Reconstituted antibody may be stored at 2 to 8 °C for up to one month. Addition of a preservative (15 mM sodium azide) may be necessary.
- For extended storage, add an equal volume of high purity glycerol, to a final concentration of 50% and BSA to a final concentration of 1% and store at -20 °C.
- During shipment, small volumes will occasionally become entrapped in the seal of the product vial. We recommend briefly centrifuging the vial in a tabletop centrifuge to dislodge any liquid in the container's cap.

Storage/Stability

Store ammonium sulfate suspension at 2-8 °C for up to one month.
For extended storage, freeze in working aliquots.
Reconstituted and diluted antiserum should be stored in aliquots at -20 °C.

Product Profile

Recommended dilutions are 1:200 to 1:1,000 for immunoblotting and immunohistochemistry.

Note: In order to obtain the best results and assay sensitivity in various techniques and preparations, we recommend determining optimal working dilutions by titration.

KAA 11/04

Sigma brand products are sold through Sigma-Aldrich, Inc.

Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.