

17119 Aeromonas Selective Supplement

Antibiotic Supplement is recommended for selective isolation of *Aeromonas* species.

Composition

(per vial, sufficient for 500 ml medium)

Ampicillin 2.5 mg.

Balanced with carbohydrates/buffer salts.

Directions

Rehydrate the content of 1 vial aseptically with 2 ml of 50% (v/v) alcohol. Mix well and aseptically add to 500 ml sterile, molten *Aeromonas* Isolation Medium Base (Cat. No. 17118). Mix well and pour into sterile petri plates.

Precautions and Disclaimer

This product is for R&D use only, not for drug, household, or other uses. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

