



RABBIT ANTI-OSTEOCALCIN, a.a. 1-22 POLYCLONAL ANTIBODY

CATALOG NUMBER:	AB1857
LOT NUMBER:	
QUANTITY:	100 μ L
SPECIFICITY:	Recognizes human osteocalcin (a.a. 1-22) mature form. Also reacts with osteocalcin (1-44) and osteocalcin (1-49). No reactivity with human bone sialoprotein and human osteonectin.
IMMUNOGEN:	Synthetic peptide corresponding to amino acids 1-22 of human Osteocalcin.
APPLICATIONS:	RIA (1:2,000) ELISA (1:5,000) Optimal working dilutions must be determined by the end user.
SPECIES REACTIVITIES:	Human. Other species have not been tested.
FORMAT:	Rabbit serum.
PRESENTATION:	Lyophilized. Reconstitute with 100 μ L of sterile distilled water.
STORAGE/HANDLING:	Maintain lyophilized material at 2-8°C. After reconstitution maintain at -20°C in undiluted aliquots for up to 6 months. Avoid repeated freeze/thaw cycles.

Important Note: *During shipment, small volumes of product will occasionally become entrapped in the seal of the product vial. For products with volumes of 200 μ L or less, we recommend gently tapping the vial on a hard surface or briefly centrifuging the vial in a tabletop centrifuge to dislodge any liquid in the container's cap.*

**FOR RESEARCH USE ONLY; NOT FOR USE IN DIAGNOSTIC
PROCEDURES. NOT FOR HUMAN OR ANIMAL CONSUMPTION**

Unless otherwise stated in our catalog or other company documentation accompanying the product(s), our products are intended for research use only and are not to be used for any other purpose, which includes but is not limited to, unauthorized commercial uses, in vitro diagnostic uses, ex vivo or in vivo therapeutic uses or any type of consumption or application to humans or animals.

©2002 - 2011: Millipore Corporation. All rights reserved. No part of these works may be reproduced in any form without permission in writing.

28820 Single Oak Drive • Temecula, CA 92590
Technical Support: T: 1-800-MILLIPORE (1-800-645-5476) • F: 1-800-437-7502
www.millipore.com

AB1857/02-SEP-2011/BKC