



MOUSE ANTI-HEPATITIS C VIRUS MONOCLONAL ANTIBODY

CATALOG NUMBER:	MAB8694	QUANTITY:	100µg
LOT NUMBER:		CONCENTRATION:	Xx mg/mL
HOST/ISOTYPE:	Ms/IgG ₁		
SPECIFICITY:	Specific to HCV NS5A		
IMMUNOGEN:	Recombinant protein to NS5A		
APPLICATIONS:	ELISA Immunofluorescence RIA Western Blot Final dilutions must be determined by end user		
SPECIES REACTIVITY:	Human		
FORMAT:	Purified immunoglobulin		
PRESENTATION:	Liquid in PBS.		
STORAGE/HANDLING:	Store at 2 to 8°C for up to 12 months from date of receipt.		

Important Note: *During shipment, small volumes of product will occasionally become entrapped in the seal of the product vial. For products with volumes of 200 µL or less, we recommend gently tapping the vial on a hard surface or briefly centrifuging the vial in a tabletop centrifuge to dislodge any liquid in the container's cap.*

FOR RESEARCH USE ONLY; NOT FOR USE IN DIAGNOSTIC
PROCEDURES. NOT FOR HUMAN OR ANIMAL CONSUMPTION

Unless otherwise stated in our catalog or other company documentation accompanying the product(s), our products are intended for research use only and are not to be used for any other purpose, which includes but is not limited to, unauthorized commercial uses, in vitro diagnostic uses, ex vivo or in vivo therapeutic uses or any type of consumption or application to humans or animals.

©2002 - 2011: Millipore Corporation. All rights reserved. No part of these works may be reproduced in any form without permission in writing.