

68075 *C. perfringens* Millichrome™ plus Supplement 2

Selective supplement for detection and enumeration of *C. perfringens*

Package information:

1 vial à 0.6g (with screw cap) is sufficient for 5L

Directions

- Add 0.12 g of *C. perfringens* Millichrome™ plus supplement 2 in 1 mL of purified water.
- Filter to sterilize at 0.45 µm
- Add 1 mL/L of the supplement 2 solution to the melted base at 45/50 °C.

Precautions and Disclaimer

This product is for R&D use only, not for drug, household, or other uses. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

