

40784 Listeria mono Selective Supplement

A selective supplement recommended by ISO Committee for the isolation of Listeria species.

Composition:

(per vial, sufficient for 500 ml medium)

Ceftazidime 10.0 mg

Nalidixic acid, sodium salt 10.0 mg

Amphotericin B 5.0 mg

Polymyxin B sulphate 38350 IU

Directions:

Rehydrate the contents of 1 vial aseptically with 2 ml of 0.2 N Sodium hydroxide, further add 8 ml of sterile distilled water. Mix well and aseptically add it to 465 ml of sterile, molten L. mono Differential Agar Base (Cat. No. 77408) along with sterile contents of Listeria mono Enrichement Supplement I (Cat. No. 03708). Mix well and pour into sterile petri plates.

Precautions and Disclaimer

This product is for R&D use only, not for drug, household, or other uses. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

