

QC RANGE SHEET
MILLIPLEX® MAP Human CVD Magnetic Bead Panel 6
Kit Catalog# HCVD6MAG-67K
QC Catalog# HCVMG6-6067-6
Lot# HCVMG6-103 and HCVMG6-203



<u>Analyte</u>	<u>QC Level</u>	<u>Expected Range</u>	<u>Units</u>
CD14	Control 1	0.15 – 0.32	ng/mL
	Control 2	1.2 – 2.5	ng/mL
DPP4	Control 1	0.79 – 1.6	ng/mL
	Control 2	6.1 – 13	ng/mL
Endostatin	Control 1	0.04 – 0.08	ng/mL
	Control 2	0.31 – 0.65	ng/mL
Progranulin	Control 1	0.39 – 0.82	ng/mL
	Control 2	3.0 – 6.3	ng/mL
LRG1	Control 1	0.36 – 0.75	ng/mL
	Control 2	3.1 – 6.4	ng/mL
MCAM	Control 1	0.41 – 0.86	ng/mL
	Control 2	3.1 – 6.4	ng/mL
ucMGP	Control 1	0.77 – 1.6	ng/mL
	Control 2	6.2 – 13	ng/mL

HCVMG6-6067-6: HCVMG6-103.203 (25 January 2021)
 EMD Millipore • 1-800-645-5476 • 1-800-645-5439 (FAX) • emdmillipore.com