


p3XFLAG-CMV-7.1 (4.7 kb)



p3XFLAG-CMV-8 (4.8 kb)



Multiple Cloning Site

(p3XFLAG-CMV-7.1 and p3XFLAG-CMV-8)

3XFLAG Peptide Sequence															
Met*	Asp	Tyr	Lys	Asp	His	Asp	Gly	Asp	Tyr	Lys	Asp	His	Asp	Ile	
ACC	ATG	GAC	TAC	AAA	GAC	CAT	GAC	GGT	GAT	TAT	AAA	GAT	CAT	GAC	ATC
TGG	TAC	CTG	ATG	TTT	CTG	GTA	CTG	CCA	CTA	ATA	TTT	CTA	GTA	CTG	TAG

3XFLAG Peptide Sequence									
Asp	Tyr	Lys	Asp	Asp	Asp	Asp	Lys		
GAT	TAC	AAG	GAT	GAC	GAT	GAC	AAG	CTT	GCG
CTA	ATG	TTC	CTA	CTG	CTA	CTG	TTC	GAA	CGC

Not I		EcoR I		Hind III	
GCG	GCC	GCG	AAT	TCA	TCG
CGC	CGG	CGC	TTA	AGT	AGC

Bgl II	EcoR V	Kpn I	Sal I	Xba I	BamH I	Sma I
ATA	GAT	CTG	ATA	TCG	GTA	CCA
TAT	CTA	GAC	TAT	AGC	CAT	GGT

*For p3XFLAG-CMV-8, the Met-preprotrypsin leader sequence (PPT LS) precedes the FLAG coding sequence.