



QUALITY CONTROL RANGES

Human PYY (Total) ELISA Catalog #EZHPYYT66K

STANDARD, LYOPHILIZED CATALOG #E8066-K

Standard	Human Total PYY Value
Lot Number: E66-20K	1827 pg/mL

The exact value of lyophilized standard is shown above.

Reconstitute with **0.5mL** distilled or deionized water.

QUALITY CONTROLS 1 AND 2, LYOPHILIZED CATALOG #E6066-K

Reconstitute each vial with **0.5mL** distilled or deionized water.

CONTROL	LOT NUMBER	EXPECTED RANGE	UNITS
QC 1	PYYT – 120	97 – 201	pg/mL
QC 2	PYYT- 220	328 – 681	pg/mL

E6066-K: PYYT-120.220 (09 August 2024)

EMD Millipore • 1-800-645-5476 • 1-800-645-5439 (FAX) • emdmillipore.com

The M logo is a trademark of Merck KGaA, Darmstadt, Germany.

© 2013 EMD Millipore Corporation, Billerica, MA 01821 USA