



## QUALITY CONTROL RANGES

### Rat FGF-23 ELISA KIT CATALOG # EZRFGF23-42K

#### STANDARD, LYOPHILIZED CATALOG # E8042-K

STANDARD	Rat FGF-23 VALUE
LOT NUMBER: 842-4K	20 ng/mL

The exact value of lyophilized standard is shown above.

Reconstitute with **0.25 mL** distilled or deionized water.

#### QUALITY CONTROLS 1 AND 2, LYOPHILIZED CATALOG # E6042-K

Reconstitute each vial with **0.25 mL** distilled or deionized water.

CONTROL	LOT NUMBER	EXPECTED RANGE	UNITS
QC 1	FGF23-104A	0.22 – 0.60	ng/mL
QC 2	FGF23-204A	2.0 – 4.1	ng/mL

EMD Millipore • (636) 441-8400 • [www.emdmillipore.com](http://www.emdmillipore.com)  
E6042-K.104A.204A (16 Aug 2016)

The M logo is a trademark of Merck KGaA, Darmstadt, Germany.  
© 2013 EMD Millipore Corporation, Billerica, MA 01821 USA.