

3050 Spruce Street Saint Louis, Missouri 63103 USA Telephone 800-325-5832 • (314) 771-5765 Fax (314) 286-7828 email: techserv@sial.com sigma-aldrich.com

ProductInformation

Synonyms: 9,10-Secoergosta-5,7,10(19),22-tetraen-

3-ol, (3β,5*Z*,7*E*,22*E*)-.

ERGOCALCIFEROL SIGMA REFERENCE STANDARD

Product No. E1007

Lot No: 79H8806

Storage Temperature: 8° - 15° or in refrigerator

Formula: C₂₈H₄₄O

Formula Weight: 396.66

Product Description

This product meets USP specifications (Current through USP 23, 10th Supplement), and is traceable to USP Reference Standard lot O


Store sealed vial in a cold place, protected from light. Allow to attain room temperature before opening ampule. Use the material promptly and discard the unused portion.


CAS No: 50-14-6

TEST	USP 23 SPECIFICATIONS	Result Sigma R.S. Lot 79H8806
Infrared absorption spectrum KBr dispersion (197K)	compares to standard	conforms to USP RS lot O
Ultraviolet absorption 10 $\mu g/mL$ in alcohol solution $\langle 197U \rangle$	compares to standard and absorptivities at 265 nm do not differ by more than 3.0%	passes
Identification C	a bright red color is produced and rapidly changes through violet and blue to green	passes
Identification D Thin layer chromatography (621)	compares to standard	passes
Melting range, °C <i>Class 1b</i> (741)	between 115° and 119°	115°-117°
Specific rotation, [α] ²⁵ 15 mg/mL in alcohol (781S)	between +103° and +106°	+104°
Reducing substances	absorbance at 525 nm is not greater than that obtained with hydroquinone, similarly treated	passes
Organic volatile impurities <i>Method V</i> (467⟩	meets requirements	passes
Assay - HPLC (621)	not less than 97.0 percent and not more than 103.0 percent	100.3%


See reverse side for IR and UV spectrum of Sigma Reference Standard.

ERGOCALCIFEROL SIGMA REFERENCE STANDARD Product No. E1007 Lot 79H8806





Wavenumbers (cm⁻¹)



7/99

Sigma brand products are sold through Sigma-Aldrich, Inc. Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.