

# **Product Information**

# Dulbecco's Modified Eagle's Medium/ **Ham's Nutrient Mixture F12**

with 3151 mg/L dextrose, with 2.5 mM L-glutamine, with 55 mg/L sodium pyruvate, without sodium bicarbonate

Description

CATALOG NO. 56498C

Studies to determine the nutritional requirements of many cells have been in progress since Eagle's first reports. The major essential nutrients were identified and work became focused on the media requirements of individual cell types. Many media designed for these purposes are now available. Among the first of these media, developed initially to study hormonal requirements of cells in culture, was a mixture of Dulbecco's Modified Eagle's Medium and Ham's Nutrient Mixture F12, known as DMEM/F12.

# **Formulation**

	• >	
Component (all components measured in mg/L)		
INORGANIC SALTS		
Calcium chloride anhydrous	116.610	
Cupric sulfate pentahydrate	0.00125	
Ferric nitrate nonahydrate	0.050	
Ferrous sulfate heptahydrate	0.417	
Magnesium chloride anhydrous	28.610	
Magnesium sulfate anhydrous	48.840	
Potassium chloride	311.800	
Sodium chloride	6999.500	
Sodium phosphate dibasic anhydrous	71.020	
Sodium phosphate monobasic monohydrate	62.500	
Zinc sulfate heptahydrate	0.4315	
VITAMINS		
Biotin	0.00365	
D-calcium pantothenate	2.240	
Choline chloride	8.980	
Cyanocobalamin	0.680	
Folic acid	2.650	
i-inositol	12.600	
Niacinamide	2.0185	
Pyridoxal HCl	2.000	
Pyridoxine HCI	0.031	
Riboflavin	0.219	
Thiamine HCI	2.170	

Component continued (all components measured in mg/L)	
AMINO ACIDS	
L-alanine	4.455
L-arginine HCI	147.500
L-asparagine monohydrate	7.500
L-aspartic acid	6.650
L-cysteine HCI monohydrate	17.560
L-cystine 2HCl	31.290
L-glutamic acid	7.350
L-glutamine	365.000
Glycine	18.750
L-histidine HCI monohydrate	31.480
L-isoleucine	54.470
L-leucine	59.050
L-lysine HCI	91.250
L-methionine	17.240
L-phenylalanine	35.480
L-proline	17.250
L-serine	26.250
L-threonine	53.450
L-tryptophan	9.020
L-tyrosine 2Na dihydrate	55.790
L-valine	52.850
OTHER	
Dextrose anhydrous	3151.000
Hypoxanthine sodium salt	2.390
Linoleic acid	0.042
DL-α-Lipoic acid	0.105
Phenol red sodium salt	8.602
Putrescine 2HCI	0.081
Sodium pyruvate	55.000
Thymidine	0.365
ADD: Sodium bicarbonate	2438.000
Grams of powder per liter:	12.001

### **Precautions**

Use aseptic technique when handling or supplementing this medium after filtration. This product is for further manufacturing use. THIS PRODUCT IS NOT INTENDED FOR HUMAN OR THERAPEUTIC USE.

# Storage

Store dry powder medium at 2 to 8 C. Do not use after expiration date. Store hydrated medium protected from light at 2 to 8 C.

### **United States**

SAFC Biosciences, Inc. 13804 W. 107th Street Lenexa, Kansas 66215

USA +1 913-469-5580 Toll free-USA 1 800-255-6032 +1 913-469-5584 E-mail

#### Europe

SAFC Biosciences Ltd. Smeaton Road, West Portway Andover, Hampshire SP10 3LF UNITED KINGDOM

Phone +44 (0)1264-333311 +44 (0)1264-332412 F-mail info-eu@sial.com

#### Asia Pacific

SAFC Biosciences Pty. Ltd. 18-20 Export Drive Brooklyn, Victoria 3025 AUSTRALIA

+61 (0)3-9362-4500 Toll free-AUS 1 800-200-404 +61 (0)3-9315-1656 E-mail info-ap@sial.com

### Indications of Deterioration

Dry powder medium should be free flowing. Do not use if powder is caked. Prepared medium should be clear of particulates and flocculent material. Do not use if liquid medium is cloudy or contains precipitate. Other evidence of deterioration may include color change or degradation of physical or performance characteristics.

# **Preparation Instructions**

- 1. Measure 80 90% of the final volume of cell culture grade water into an appropriate size mixing vessel. Water temperature should be 15 to 30 C. Do not heat water.
- Add the dry powder medium to the water slowly. Rinse the original package with a small amount of cell culture grade water to remove all traces of powder and add to the solution. Mix until completely dissolved.
- 3. For each liter of DMEM/F12 being prepared, add 2.44 g/L of sodium bicarbonate (Catalog No. 90421C) or 32.5 mL/L of sodium bicarbonate solution 7.5% (Catalog No. 59221C) solution. Mix until completely dissolved.
- 4. While stirring the solution in Step 3, adjust the pH to 6.9 7.1 using NaOH 1N (Catalog No. 59223C) or HCl 1N. The pH of bicarbonate buffered solutions usually rises 0.1 0.2 units during filtration.
- 5. Add cell culture grade water to the solution in Step 4 to bring it to the final volume. Keep the vessel closed until the solution is filtered to avoid fluctuations in pH.
- 6. Sterilize the solution using a membrane filter with a pore size of 0.2  $\mu$ m or less. A peristaltic pump or an inert gas such as nitrogen can be used to provide positive pressure at 3 15 psi. Do not use CO<sub>2</sub> gas.
- 7. Sterile solutions should be dispensed aseptically into sterile containers. Store protected from light at 2 to 8 C.
- 8. Supplements, such as antibiotics, can be added to the sterilized solution using aseptic technique. Storage conditions and shelf life of the supplemented product may be affected by the nature of the supplements. Sterile serum should not be refiltered before or after being added to sterile medium because growth promoting capacity may be reduced upon refiltration.

**NOTE:** Dry powder medium is extremely hygroscopic and must be protected from atmospheric moisture. We recommend that the entire contents of each package be used immediately after opening.

## Characteristics

**Appearance** 

Off-white free-flowing powder **Bioburden** 

≤ 100 CFU/100 mL

**Endotoxin** 

≤ 1.0 EU/mL

Osmolality (as supplied)

240 - 280 mOsm/kg H<sub>2</sub>O

Osmolality (with NaHCO<sub>3</sub>)

285 - 325 mOsm/kg H<sub>2</sub>O

pH (as supplied)

6.3 - 6.7

pH (with NaHCO<sub>3</sub>)

7.5 - 7.0

#### Warranty, Limitation of Remedies

SAFC Biosciences warrants to the purchaser for a period of one year from date of delivery that this product conforms to its specifications. Other terms and conditions of this warranty are contained in SAFC Biosciences' written warranty, as copy of which is available upon request. ALL OTHER WARRANTIES, EXPRESSED DOR IMPLIED, INCLUDING THE IMPLIED WARRANTY OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE, ARE EXCLUDED. In no case will SAFC Biosciences be liable for any special, incidental, or consequential damages arising out of this product or the use of this product by the customer or any third party based upon breach of warranty, breach of contract, negligence, strict tort, or any other legal theory. SAFC Biosciences expressly disclaims any warranty against claims by any third party by way of infringement or the like. THIS PRODUCT IS INTENDED FOR PURPOSES DESCRIBED ONLY AND IS NOT INTENDED FOR ANY HUMAN OR THERAPEUTIC USE.

Additional Terms and Conditions are contained in the product Catalog, a copy of which is available upon request

© 2006 SAFC Biosciences, Inc.

Issued December 2006 P56498 1203 1005 0406 0606 0906

**United States** 

SAFC Biosciences, Inc. 13804 W. 107th Street Lenexa, Kansas 66215

USA Phone Toll free-USA

+1 913-469-5580 A 1 800-255-6032 +1 913-469-5584 Europe

SAFC Biosciences Ltd. Smeaton Road, West Portway Andover, Hampshire SP10 3LF UNITED KINGDOM

UNITED KINGDOM
Phone +44 (0)1264-333311
Fax +44 (0)1264-332412
E-mail info-eu@sial.com

Asia Pacific

SAFC Biosciences Pty. Ltd. 18-20 Export Drive Brooklyn, Victoria 3025 AUSTRALIA

AUSTRALIA
Phone +61 (0)3-9362-4500
Toll free-AUS 1 800-200-404
Fax +61 (0)3-9315-1656