



QUALITY CONTROL RANGES
Human Adipocyte Magnetic Bead Panel
Catalog # HADCYMAG-61K
QC Catalog #: HADCYT-6061

Control 1 Lot#: HADCY-110 / Control 2 Lot#: HADCY-210

<u>Analyte</u>	<u>QC Level</u>	<u>Expected Range</u>	<u>Units</u>
NGF	Control 1	27 - 55	pg/mL
	Control 2	125 - 260	pg/mL
IL-6	Control 1	72 - 150	pg/mL
	Control 2	394 - 819	pg/mL
LEPTIN	Control 1	869 - 1804	pg/mL
	Control 2	4238 - 8802	pg/mL
IL-8	Control 1	20 - 42	pg/mL
	Control 2	103 - 215	pg/mL
HGF	Control 1	527 - 1095	pg/mL
	Control 2	2517 - 5228	pg/mL
ADPN	Control 1	1761 - 3657	pg/mL
	Control 2	9639 - 20019	pg/mL
MCP-1	Control 1	87 - 180	pg/mL
	Control 2	389 - 807	pg/mL
TNF- α	Control 1	32 - 66	pg/mL
	Control 2	153 - 317	pg/mL
RESISTIN	Control 1	378 - 785	pg/mL
	Control 2	2084 - 4328	pg/mL
IL-1B	Control 1	41 - 85	pg/mL
	Control 2	208 - 432	pg/mL
PAI-1(total)	Control 1	776 - 1612	pg/mL
	Control 2	3899 - 8098	pg/mL

HADCYT-6061: HADCY-110.210 (12 February 2025)

EMD Millipore • 1-800-645-5476 • 1-800-645-5439 (FAX) • www.emdmillipore.com

The M logo is a trademark of Merck KGaA, Darmstadt, Germany.

© 2013 EMD Millipore Corporation, Billerica, MA 01821 USA