



QUALITY CONTROL RANGES
Human Adipocyte Magnetic Bead Panel
Catalog # HADCYMAG-61K
QC Catalog #: HADCYT-6061

Control 1 Lot#: HADCY-109 / Control 2 Lot#: HADCY-209

<u>Analyte</u>	<u>QC Level</u>	<u>Expected Range</u>	<u>Units</u>
NGF	Control 1	29 - 60	pg/mL
	Control 2	132 - 274	pg/mL
IL-6	Control 1	79 - 164	pg/mL
	Control 2	385 - 799	pg/mL
LEPTIN	Control 1	916 - 1902	pg/mL
	Control 2	4316 - 8963	pg/mL
IL-8	Control 1	21 - 44	pg/mL
	Control 2	105 - 218	pg/mL
HGF	Control 1	562 - 1168	pg/mL
	Control 2	2466 - 5122	pg/mL
ADPN	Control 1	1973 - 4098	pg/mL
	Control 2	10563 - 21938	pg/mL
MCP-1	Control 1	88 - 182	pg/mL
	Control 2	391 - 811	pg/mL
TNF- α	Control 1	32 - 67	pg/mL
	Control 2	161 - 334	pg/mL
RESISTIN	Control 1	560 - 1163	pg/mL
	Control 2	2648 - 5499	pg/mL
IL-1B	Control 1	43 - 89	pg/mL
	Control 2	213 - 443	pg/mL
PAI-1(total)	Control 1	818 - 1699	pg/mL
	Control 2	3906 - 8112	pg/mL

HADCYT-6061: HADCY-109.209 (11 December 2023)

EMD Millipore • 1-800-645-5476 • 1-800-645-5439 (FAX) • www.emdmillipore.com

The M logo is a trademark of Merck KGaA, Darmstadt, Germany.

© 2013 EMD Millipore Corporation, Billerica, MA 01821 USA