


## CHEMISCREEN™ MEMBRANE PREPARATION RECOMBINANT HUMAN GPR103/QRFP NEUROPEPTIDE RECEPTOR


<b>CATALOG NUMBER:</b>	HTS189M	<b>QUANTITY:</b>	200 units
<b>LOT NUMBER:</b>	R0708E0057	<b>VOLUME/CONCENTRATION:</b>	2 mL, 1 mg/mL

**BACKGROUND:** GPR103 is a Gq-coupled receptor that is located in the brain and adrenal gland, and is activated by a 43-amino acid RFamide peptide, QRFP43, and a 26 amino acid fragment, QRFP26/26RFa/Peptide 518 (Jiang *et al.*, 2003; Fukusumi *et al.*, 2003). QRFP43 is the active form in rat brain, and intracerebral administration of QRFP43 induces feeding and metabolic rate (Takayasu *et al.*, 2006; Moriya *et al.*, 2006). Millipore's GPR103 membrane preparations are crude membrane preparations made from our proprietary stable recombinant cell lines to ensure high-level of GPCR surface expression; thus, they are ideal HTS tools for screening of antagonists of GPR103. The membrane preparations exhibit a K<sub>d</sub> of 0.63nM for [<sup>125</sup>I]-QRFP43. With 10 µg/well GPR103 Membrane Prep and 0.35 nM [<sup>125</sup>I]- QRFP43, a greater than 6-fold signal-to-background ratio was obtained.

**APPLICATIONS:** Radioligand binding assay



**Figure 1. Saturation binding for GPR103.** 10 µg/well GPR103 Membrane Preparation was incubated with increasing amount of [<sup>125</sup>I]QRFP43 in the absence (total binding, TB) or presence (nonspecific binding, NSB) of 500-fold excess unlabeled 26Ffa human Hypothalamic Peptide. Specific binding (SB) was determined by subtracting NSB from TB.



**Figure 2. Competition binding for GPR103.** 10  $\mu$ g/well GPR103 Membrane Preparation or Wild-Type Chem-1 membrane preparation (WT; Chemicon Catalog # HTS000MC1) was incubated with 0.35 nM [ $^{125}$ I]-QRFP43 and increasing concentrations of unlabeled 25Rfa, and more than 6- fold signal:background was obtained.

**Table 1.** Signal:background and specific binding values obtained in a competition binding assay with varying amounts of GPR103 membrane prep.

	10 $\mu$ g/well
Signal:background	10
Specific binding (cpm)	4217.1

SPECIFICATIONS: 1 unit = 10  $\mu$ g membrane preparation  
Bmax: 0.81 pmol/mg  
K<sub>d</sub>: 0.63 nM

Species: Human GPR103 (Accession Number: NM\_198179)

HOST CELLS: Chem-1, an adherent mammalian cell line without any endogenous GPR103 expression.

RECOMMENDED ASSAY CONDITIONS: Membranes are mixed with radioactive ligand and unlabeled competitor (see Figures 1 and 2 for concentrations tested) in binding buffer in a nonbinding 96-well plate, and incubated for 1-2 h. Prior to filtration, an FC 96-well harvest plate (Millipore cat. # MAHF C1H) is coated with 0.33% polyethylenimine for 30 min, then washed with 50mM HEPES, pH 7.4, 0.5% BSA. Binding reaction is transferred to the filter plate, and washed 3 times (1 mL per well per wash) with Wash Buffer. The plate is dried and counted.

Binding buffer: 50 mM Hepes, pH 7.4, 5 mM MgCl<sub>2</sub>, 1 mM CaCl<sub>2</sub>, 0.2% BSA, filtered and stored at 4°C

Radioligand: [ $^{125}$ I]QRFP43 (Perkin Elmer # NEX408)

Wash Buffer: 50 mM Hepes, pH 7.4, 500mM NaCl, 0.1% BSA, filtered and stored at 4°C.

One package contains enough membranes for at least 200 assays (units), where an unit is the amount of membrane that will yield greater than 6-fold signal:background with  $^{125}\text{I}$ -labeled QRFP43 at 0.35 nM.

**PRESENTATION:** Liquid in packaging buffer: 50 mM Tris pH 7.4, 10% glycerol and 1% BSA with no preservatives.  
Packaging method: Membrane protein was adjusted to the indicated concentration in packaging buffer, rapidly frozen, and stored at  $-80^{\circ}\text{C}$ .

**STORAGE/HANDLING:** Maintain frozen at  $-70^{\circ}\text{C}$  for up to 2 years. Do not freeze and thaw.

**REFERENCES:** Fukusumi S *et al.* (2003) A new peptidic ligand and its receptor regulating adrenal function in rats. *J. Biol. Chem.* 278: 46387-46395.

Jiang Y. *et al.* (2003) Identification and characterization of a novel RF-amide peptide ligand for orphan G-protein-coupled receptor SP9155. *J. Biol. Chem.* 278: 27652-27657.

Moriya R *et al.* (2006) RFamide peptide QRFP43 causes obesity with hyperphagia and reduced thermogenesis in mice. *Endocrinology* 147: 2916-2922.

Takayasu S *et al.* (2006) A neuropeptide ligand of the G protein-coupled receptor GPR103 regulates feeding, behavioral arousal, and blood pressure in mice. *Proc. Natl. Acad. Sci USA* 103: 7438-7443.

*For research use only; not for use as a diagnostic.*

Unless otherwise stated in our catalog or other company documentation accompanying the product(s), our products are intended for research use only and are not to be used for any other purpose, which includes but is not limited to, unauthorized commercial uses, in vitro diagnostic uses, ex vivo or in vivo therapeutic uses or any type of consumption or application to humans or animals.

©2007: Millipore Corporation. All rights reserved. No part of these works may be reproduced in any form without permission in writing.