

Product Information

Monoclonal Anti-FLAG®-Peroxidase produced in rat

Clone 6F7, purified immunoglobulin

SAB4200119

Product Description

Epitope tags provide a method to localize gene products in a variety of cell types, study the topology of proteins and protein complexes, identify associated proteins, and characterize newly identified, low abundance, or poorly immunogenic proteins when protein-specific antibodies are not available. The FLAG® peptide sequence, known also as DYKDDDDK, is one of the most widely used protein tags in recombinant protein expression and purification.¹ Protein tagging with the FLAG® tag may be done at the N-terminus, the N-terminus preceded by a methionine residue, the C-terminus, or at internal positions of the target protein. The small size of the FLAG® tag or sequence and its high hydrophilicity tend to decrease the possibility of interference with the protein expression, proteolytic maturation, antigenicity, and function.

Monoclonal Anti-FLAG®-Peroxidase is a purified immunoglobulin fraction of monoclonal Anti-FLAG® (rat IgG1 isotype) isolated from culture supernatant of 6F7 hybridoma cells grown in a bioreactor, conjugated to horseradish peroxidase (HRP). The hybridoma 6F7 was produced by the fusion of mouse myeloma cells and splenocytes from rat which was immunized with the FLAG peptide.

Monoclonal Anti-FLAG®-Peroxidase, recognizes N-terminal, C-terminal and internal FLAG®-tagged fusion proteins. This product is especially recommended for identifying C-terminal FLAG®-tagged fusion proteins. The product can be used for immunoblotting. Several publications cite use of product SAB4200119 in their protocols.²⁻⁵

Reagent

This product is supplied as a solution in 0.01 M phosphate buffered saline, pH 7.4, containing 0.01% merthiolate as a preservative.

Antibody concentration

2-4 mg/mL (exact value on Certificate of Analysis for particular lot)

Molar ratio Ab/E

0.6-1.5 (exact value on Certificate of Analysis for particular lot)

Storage/Stability

- For continuous use, store at 2-8 °C for up to one month.
- For extended storage, freeze at -20 °C in working aliquots.
- Repeated freezing and thawing, or storage in "frost-free" freezers, is not recommended.
- If slight turbidity occurs upon prolonged storage, clarify the solution by centrifugation before use.
- Working dilution samples should be discarded if not used within 12 hours.

Product Profile

Immunoblotting

A working dilution of 1:1000-1:2000 is recommended, using extracts of transfected cells that express C-terminal FLAG®-tagged fusion protein.

Note: To obtain best results and assay sensitivity with various techniques and preparations, we recommend determining optimal working dilutions by titration.

Procedure

Procedure for Immunoblotting

1. Separate FLAG®-tagged fusion proteins from sample lysates using a standard SDS-PAGE protocol. Load 2.5-20 µg of total lysate protein per lane.
Note: The amount of extract to be loaded per slab or lane depends on the level of protein expression and may vary between experiments.
2. Transfer proteins from the gel to a nitrocellulose membrane.
3. Block the membrane using a solution of 5% non-fat dry milk in Dulbecco's Phosphate Buffered Saline (DPBS, Cat. No. D8537) at room temperature for 1 hour.
4. Wash the membrane three times for 5 minutes each in PBS containing 0.05% TWEEN® 20 (PBS-T, Cat. No. P3563) at room temperature.
5. Incubate the membrane with Monoclonal Anti-FLAG®-Peroxidase using an optimized concentration in PBS containing 0.5% non-fat dry milk, at room temperature with agitation for 2 hours.
6. Wash the membrane three times for 5 minutes each in PBS containing 0.05% TWEEN® 20 at room temperature.
7. Treat the membrane with a peroxidase substrate.

References

1. Terpe, K., *Appl. Microbiol. Biotechnol.*, **60(5)**, 523-533 (2003).
2. Libby, E.A. *et al.*, *PLoS Genet.*, **11(6)**, e1005275 (2015).
3. Valsecchi, C.I.K. *et al.*, *Nat. Commun.*, **9(1)**, 3626 (2018).
4. Lee, D.H. *et al.*, *Nat. Commun.*, **11(1)**, 1838 (2020).
5. Valsecchi, C.I.K. *et al.*, *Nature*, **589(7840)**, 137-142 (2021).

Notice

We provide information and advice to our customers on application technologies and regulatory matters to the best of our knowledge and ability, but without obligation or liability. Existing laws and regulations are to be observed in all cases by our customers. This also applies in respect to any rights of third parties. Our information and advice do not relieve our customers of their own responsibility for checking the suitability of our products for the envisaged purpose.

The information in this document is subject to change without notice and should not be construed as a commitment by the manufacturing or selling entity, or an affiliate. We assume no responsibility for any errors that may appear in this document.

Technical Assistance

Visit the tech service page at SigmaAldrich.com/techservice.

Standard Warranty

The applicable warranty for the products listed in this publication may be found at SigmaAldrich.com/terms.

Contact Information

For the location of the office nearest you, go to SigmaAldrich.com/offices.

The life science business of Merck KGaA, Darmstadt, Germany operates as MilliporeSigma in the U.S. and Canada.

MilliporeSigma, and Sigma-Aldrich are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.
© 2021 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved.

SAB4200119dat Rev 09/21 ST,TD,KAA,PHC,GCY

For research use only. Not for use in diagnostic procedures.