3050 Spruce Street, St. Louis, MO 63103 USA Tel: (800) 521-8956 (314) 771-5765 Fax: (800) 325-5052 (314) 771-5757 email: techservice@sial.com sigma-aldrich.com

# **Product Information**

## Ganciclovir

Catalog Number **G2536** Storage Temperature 2-8 °C

CAS RN: 82410-32-0

Synonyms: 9-(1,3-Dihydroxy-2-propoxymethyl)guanine; 2'-Nor-2'-deoxyguanosine; 2'-NDG; BIOLF-62; DHPG; BW-759U

# **Product Description**

Molecular formula: C<sub>9</sub>H<sub>13</sub>N<sub>5</sub>O<sub>4</sub> Molecular weight: 255.23

Ganciclovir, a nucleoside analog, can be enzymatically phosphorylated to an active triphosphate analog by herpes simplex virus thymidine kinase (HSV-TK). The phosphorylated product causes inhibition of viral DNA polymerase. Ganciclovir has also been used to study sensitivity of different viruses to antiviral treatment.<sup>1,2</sup>

This enzyme-prodrug reaction is being used to study gene therapy in cancer research. Styression of a suicide gene encoding a viral enzyme such as TK results in the phosphorylation of the non-toxic prodrug, ganciclovir. The active, phosphorylated analog is incorporated into the DNA of replicating bystander eukaryotic cells causing death of the malignant dividing cells.

Ganciclovir was found to be more cytotoxic than other HSV-TK substrates, acyclovir and araT (1- $\beta$ -D- arabinofuranosylthymine). The mechanism presumably proceeds via irreversible cell cycle arrest at the G2-M checkpoint. The involvement of gap junctions with the bystander effect observed with the HSV-TK gene and ganciclovir was reported.

The HSV-TK gene and ganciclovir are used in molecular biology for negative selection against random recombination events when homologous recombination of a gene of interest is required. Negative selection with the HSV-TK gene and ganciclovir was recently used for measuring loss of telomeres in the study of chromosome healing. Recently it has been widely used as a potent treatment against Cytomegalovirus infection following tissue transplantation. 11,12

#### **Precautions and Disclaimer**

This product is for R&D use only, not for drug, household, or other uses. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

# **Preparation Instructions**

Soluble at 10 mg/ml in 0.1 N HCl.

#### Storage/Stability

Store desiccated at 2-8  $^{\circ}$ C. Under these conditions the product is stable for 3 years.

### References

- Oon, C.J., et al., Hepatitis B virus variants with lamivudine-related mutations in the DNA polymerase and the 'a' epitope of the surface antigen are sensitive to ganciclovir., Antiviral Res., 41, 113-8 (1999).
- Cannon, J.S., et al., Human herpes virus 8-encoded thymidine kinase and phosphotransferase homologues confer sensitivity to ganciclovir., J. Virol., 73, 4786-93 (1999).
- Maatta, A.M., et al., Non-small cell lung cancer as a target disease for herpes simplex type 1 thymidine kinase-ganciclovir gene therapy., Int J. Oncol., 24, 943-9 (2004).
- 4. Wang, J., et al., Herpes simplex virus thymidine kinase and ganciclovir suicide gene therapy for human pancreatic cancer., World J. Gastroenterol., **10**, 400-3 (2004).

- Rubsam, L.Z., et al., Cytotoxicity and accumulation of ganciclovir triphosphate in bystander cells cocultured with herpes simplex virus type 1 thymidine kinase-expressing human glioblastoma cells., Cancer Res., 59, 669-75 (1999).
- Rubsam, L.Z., et al., Superior cytotoxicity with ganciclovir compared with acyclovir and 1-β-Darabinofuranosylthymine in herpes simplex virusthymidine kinase-expressing cells: a novel paradigm for cell killing., Cancer Res., 58, 3873-82 (1998).
- Halloran, P.J., and Fenton, R.G., Irreversible G2-M arrest and cytoskeletal reorganization induced by cytotoxic nucleoside analogues., Cancer Res., 58, 3855-65 (1998).
- 8. Yamasaki, H., et al., Role of connexin (gap junction) genes in cell growth control and carcinogenesis., C. R. Acad. Sci. III, **322**, 151-9 (1999).

- Introduction of DNA into Mammalian Cells, in Current Protocols in Molecular Biology, Ausubel, F.M., et al., Eds., John Wiley & Sons, Inc. (New York, NY: 1999), Volume I, Supplement 45, 9.15.1-9.15.10.
- 10. Sprung, C.N., et al., Chromosome healing in mouse embryonic stem cells., Proc. Natl. Acad. Sci. USA, **96**, 6781-6 (1999).
- Said, T., et al., Cytomegalovirus prophylaxis with ganciclovir in kidney transplant recipients receiving induction antilymphocyte antibodies., Transplant Proc., 36, 1847-9 (2004).
- Zamora, M.R., et al., Following universal prophylaxis with intravenous ganciclovir and cytomegalovirus immune globulin, valganciclovir is safe and effective for prevention of CMV infection following lung transplantation., Am. J. Transplant., 4, 1635-42 (2004).

NDH,PHC 11/10-1