

68022 *Y. enterocolitica* Millichrome™ plus Supplement

Selective supplement for *Yersinia enterocolitica*

Composition per vial:
(1 vial is sufficient for 5 L of final media)
Selective mix 0.5g

Directions

- Prepare a stock solution of the supplement: Add 100 mg to 1 mL of purified water.
- Swirl well until complete dissolution. Filter sterilize at 0.45 µm.
- Warning 2: This supplement stock solution should be used immediately after preparation, or can be stored at -20 °C and used within 15 days.
- Add 1 mL of the prepared supplement solution to 1 L of the prepared *Y. enterocolitica* Millichrome™ plus Agar base cooled at 45-50 °C +/- 2 °C.
- Swirl gently to homogenize.
- Pour into sterile Petri dishes.

Precautions and Disclaimer

This product is for R&D use only, not for drug, household, or other uses. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

