



QUALITY CONTROL RANGES

HIGH SENSITIVITY GLUCAGON-LIKE PEPTIDE (GLP-1) ACTIVE CHEMILUMINESCENT ELISA CATALOG # EZGLPHS-35K

STANDARD, LYOPHILIZED CATALOG # E8035-HS

Reconstitute standard with **0.5mL** distilled or deionized water.

STANDARD	Active GLP1 Value
LOT NUMBER: 835HS-9K	89 pM

QUALITY CONTROLS 1 AND 2, LYOPHILIZED CATALOG # E6035-HS

Reconstitute each vial with **0.5mL** distilled or deionized water.

CONTROL	LOT NUMBER	EXPECTED RANGE	UNITS
QC 1	GLPHS-109	2.0 – 4.2	pM
QC 2	GLPHS-209	13 - 26	pM

E6035-HS: GLPHS-109.209 (18 May 2022)

EMD Millipore • 1-800-645-5476 • 1-800-645-5439 (FAX) • emdmillipore.com

The M logo is a trademark of Merck KGaA, Darmstadt, Germany.

© 2013 EMD Millipore Corporation, Billerica, MA 01821 USA