



Product Information

D-(+)-Glucose Plant Cell Culture Tested

Product Number **G 7520**
Store at Room Temperature

Product Description

Molecular Formula: $C_6H_{12}O_6$
Molecular Weight: 180.2
CAS Number: 50-99-7
Melting Point : 146 °C (α -D-glucose), 150 °C (β -D-glucose)¹
pH: 5.9 (0.5 M aqueous)²

Density of solutions in water at 17.5 °C:²

Concentration (w/v%)	5	10	20	30	40
Density (g/ml)	1.019	1.038	1.076	1.113	1.149

This product is plant cell culture tested (15 mg/ml) and is appropriate for use in plant cell culture experiments.

Glucose is a main source of energy for living organisms. It occurs naturally in the free state in fruits and other parts of plants. Glucose is combined into glucosides, disaccharides, oligosaccharides, the polysaccharides (cellulose and starch), and glycogen.

Glucose is a mixture of α - and β -anomers, primarily the α -anomer. The optical rotation of the α -anomer is +112.2° (c = 100 mg/ml, 20 °C) and the β -anomer is +18.7° (c = 100 mg/ml, 20 °C). When D-glucose is dissolved in water, the optical rotation gradually changes (mutarotates) with time and approaches a final equilibrium value of +52.7° (c = 100 mg/ml, 20 °C) due to the formation of an equilibrium mixture consisting of approximately one-third α - and two-thirds β -D-glucose.²

Normal human blood contains 0.08 - 0.1% glucose.² Small amounts of glucose (also hydrogen peroxide or glucose oxidase) can be measured using luminol as a substrate with horseradish peroxidase.³

Precautions and Disclaimer

For Laboratory Use Only. Not for drug, household or other uses.

Preparation Instructions

This product is soluble in water (133 mg/ml), yielding a clear, colorless solution. Other solubility values in water have been reported as 909 mg/ml (25 °C) and 5.55 g/ml (90 °C).²

References

1. Biochemistry, 2nd ed., Lehninger, A. L., ed., Worth Publishers, Inc. (New York, NY: 1975), p. 253.
2. The Merck Index, 13th Ed., Entry# 4472.
3. Puget, K., and Michelson, A. M., Microestimation of glucose and glucose oxidase. *Biochimie*, **58**, 757-758 (1976).
4. Martindale The Extra Pharmacopoeia, 29th ed., Reynolds, J. E. F., ed., The Pharmaceutical Press (London, England: 1989), p. 1265.

MES/GCY/AJH 5/03

Sigma brand products are sold through Sigma-Aldrich, Inc.

Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.