

55286 VRE ChromoSelect Agar Supplement

VRE *ChromoSelect* Agar Supplement is recommended for selective isolation of Vancomycin Resistant Enterococci (VRE).

Composition
(per vial, sufficient for 500 ml medium)
Vancomycin 4mg
Fluconazole 5mg

Directions
Rehydrate the contents of 1 vial aseptically with 5 ml sterile distilled water. Mix well and aseptically add to 500ml of sterile, molten, cooled (45-50°C) VRE *ChromoSelect* Agar Base (Cat. No. 79599). Mix well and pour into sterile Petri plates.

Precautions and Disclaimer

This product is for R&D use only, not for drug, household, or other uses. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

