

91986 PALCAM *Listeria* Selective Supplement according to Van Netten et al.

An antimicrobial supplement recommended for the selective isolation and identification of *Listeria monocytogenes*.

Composition

(per vial, sufficient for 500 ml medium)

Ceftazidime 10.0 mg

Polymyxin B 50000 IU

Acriflavin Hydrochloride 2.5 mg

Directions

Aseptically add 5 ml of sterile distilled water to the vial. Shake to dissolve and aseptically add to 500 ml sterile, molten PALCAM *Listeria* Selective Agar (Cat. No. 75977) / PALCAM *Listeria* Selective Enrichment Broth, Vegitone (Cat. No. 59859) / PALCAM *Listeria* Selective Agar, Vegitone (Cat. No. 15776).

Mix well and dispense as desired.

Precautions and Disclaimer

This product is for R&D use only, not for drug, household, or other uses. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

