



**Human Circulating Cancer Biomarker
MILLIPLEX[®] MAP Magnetic Bead Panel 3
Kit Catalog # HCCBP3MAG-58K
Control Catalog # HCC3-6058-3
Lot# HCC3-111 and HCC3-211**

<u>Analyte</u>	<u>QC Level</u>	<u>Expected Range</u>	<u>Units</u>
Galectin	Control 1	0.016 - 0.033	ng/mL
	Control 2	0.39 - 0.80	ng/mL
Cathepsin D	Control 1	1.3 - 2.7	ng/mL
	Control 2	21 - 44	ng/mL
FAP	Control 1	0.32 - 0.66	ng/mL
	Control 2	6.4 - 13	ng/mL
MIA	Control 1	0.47 - 1.0	ng/mL
	Control 2	9.2 - 19	ng/mL
MPO	Control 1	0.13 - 0.26	ng/mL
	Control 2	2.5 - 5.1	ng/mL
SHBG	Control 1	0.14 - 0.29	nM
	Control 2	2.9 - 6.1	nM
IGFBP3	Control 1	0.37 - 0.77	ng/mL
	Control 2	8.8 - 18	ng/mL
Ferritin	Control 1	0.15 - 0.31	ng/mL
	Control 2	3.0 - 6.3	ng/mL

HCC3-6058-3: HCC3-111.211 (20 June 2024)

EMD Millipore • 1-800-645-5476 • 1-800-645-5439 (FAX) • emdmillipore.com

The M logo is a trademark of Merck KGaA, Darmstadt, Germany.

© 2013 EMD Millipore Corporation, Billerica, MA 01821 USA