



QUALITY CONTROL RANGES

MILLIPLEX[®] MAP Human Liver Injury Magnetic Bead Panel Kit Catalog # HLINJMAG-75K

Control Catalog # HLINJ-6075 Lot# HLI-106 and HLI-206

<u>Analyte</u>	<u>QC Level</u>	<u>Expected Range</u>	<u>Units</u>
ARG1	Control I	454 – 942	pg/mL
	Control II	3220 – 6800	pg/mL
MDH1	Control I	420 – 872	pg/mL
	Control II	3107 – 6453	pg/mL
GSTa	Control I	95 – 198	pg/mL
	Control II	688 – 1430	pg/mL
5-NT	Control I	248 – 515	pg/mL
	Control II	1865 – 3874	pg/mL
SDH	Control I	3108 – 7429	pg/mL
	Control II	22438 – 55575	pg/mL

HLINJ-6075: HLI-106.206 (31 December 2024)

EMD Millipore • 1-800-645-5476 • 1-800-645-5439 (FAX) • emdmillipore.com

The M logo is a trademark of Merck KGaA, Darmstadt, Germany.

© 2013 EMD Millipore Corporation, Billerica, MA 01821 USA