



Quality Control Ranges
Milliplex PLEXpedition Magnetic Bead Screening Panel
Catalog# HPLX1-9010-02
Control Catalog # HPLX1-6010-02 Lot # PLX12-101 and PLX12-201

<u>Bead</u>	<u>QC Level</u>	<u>Expected Range</u>	<u>Units</u>	<u>Bead</u>	<u>QC Level</u>	<u>Expected Range</u>	<u>Units</u>
HIL1A	Control 1	14 - 29	pg/mL	HSCF	Control 1	103 - 214	pg/mL
	Control 2	54 - 112	pg/mL		Control 2	468 - 973	pg/mL
HIL7	Control 1	9 - 18	pg/mL	HSFAS	Control1	428 - 889	pg/ml
	Control 2	36 - 74	pg/mL		Control 2	1875 - 3893	pg/ml
HIL11	Control 1	119 - 248	pg/mL	HTRAIL	Control1	105 - 218	pg/ml
	Control 2	459 - 953	pg/mL		Control 2	405 - 842	pg/ml
HIL12P40	Control 1	336 - 697	pg/mL	HIL3	Control1	28 - 58	pg/ml
	Control 2	1320 - 2743	pg/mL		Control 2	111 - 231	pg/ml
HIL12P70	Control 1	53 - 110	pg/mL	HIL9	Control1	134 - 278	pg/ml
	Control 2	219 - 455	pg/mL		Control 2	517 - 1074	pg/ml
HIL18	Control 1	6 - 12	pg/mL	HIL17F	Control1	28 - 59	pg/ml
	Control 2	24 - 50	pg/mL		Control 2	109 - 227	pg/ml
HIL22	Control 1	248 - 516	pg/mL	HMDC	Control1	23 - 48	pg/ml
	Control 2	986 - 2048	pg/mL		Control 2	89 - 186	pg/ml
HFLT3L	Control 1	20 - 42	pg/mL	HTNFB	Control1	56 - 116	pg/ml
	Control 2	77 - 161	pg/mL		Control 2	226 - 470	pg/ml
HPRFRN	Control 1	500 - 1038	pg/mL	HTSLP	Control1	3 -6	pg/ml
	Control 2	2152 - 4469	pg/mL		Control 2	11 -23	pg/ml
HCD40L	Control 1	155 - 322	pg/mL				
	Control 2	614 - 1276	pg/mL				
HGRNZMA	Control 1	29 - 61	pg/mL				
	Control 2	115 - 239	pg/mL				

Note: The Quality Control Ranges are generated with overnight assay format using serum matrix provided in the kit.
HPLX1-6010-02: PLX12-101.201 (09 April 2025) EMD Millipore • (800) 645-5476 • (800) 645-5439 (FAX) • www.millipore.com

The M logo is a trademark of Merck KGaA, Darmstadt, Germany. © 2013 EMD Millipore Corporation, Billerica, MA 01821 USA