

52099 KPC *ChromoSelect* Agar Supplement

KPC *ChromoSelect* Agar Supplement is recommended for the detection of gram negative bacteria with a reduced susceptibility to carbapenem agents.

Composition (per vial sufficient for 500 ml medium)
Selective Mix 0.200g

Directions:

Rehydrate the contents of 1 vial aseptically with 5 ml sterile distilled water. Mix well and aseptically add to 500ml of sterile, molten, cooled (45-50°C) KPC *ChromoSelect* Agar Base (Cat. No. 44022).

Mix well and pour into sterile Petri plates.

Precautions and Disclaimer

This product is for R&D use only, not for drug, household, or other uses. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

