

## 59746 Atto Rho14

### Application


Atto Rho14 is a new rhodamine featuring, a functionality for coupling to biomolecules such as DNA, RNA or proteins. The label shows strong absorption, and extraordinarily high fluorescence quantum yield. In fact Atto Rho14 is the brightest label available in this wavelength range. In addition the dye exhibits an exceptionally high photostability.

After coupling to a substrate Atto Rho14 carries a net electrical charge of +1.


Absorption and fluorescence are pH-independent in the range of pH 2 to 11, used in typical applications. The dye is moderately hydrophilic.

### Product Description

MW	884 g/mol
$\lambda_{\text{abs}}$	626 nm
$\epsilon_{\text{max}}$	$1.4 \times 10^5 \text{ M}^{-1} \text{ cm}^{-1}$
$\lambda_{\text{fl}}$	646 nm
$\eta_{\text{fl}}$	80 %
$\tau_{\text{fl}}$	3.7 ns
CF <sub>260</sub>	0.26
CF <sub>280</sub>	0.47



### Optical data of the carboxy derivative (in aqueous solution)



**Storage:** store at  $\leq -20^\circ\text{C}$ . Protect from long-term exposure to moisture and light.

**Precautions and Disclaimer** This product is for R&D use only, not for drug, household, or other uses. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

The vibrant M and Sigma-Aldrich are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. Detailed information on trademarks is available via publicly accessible resources.  
© 2018 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved.

The life science business of Merck KGaA, Darmstadt, Germany operates as MilliporeSigma in the US and Canada.

