



MOUSE ANTI-CONNEXIN 35/36 MONOCLONAL ANTIBODY

CATALOG NUMBER:	MAB3045
LOT NUMBER:	
QUANTITY:	100 µg
CONCENTRATION:	1 mg/mL
SPECIFICITY:	Recognizes perch Connexin 35. Does not cross-react with Connexin 34.7, a very closely related protein. Also recognizes rabbit and macaque Connexin 36, the gene homologue of Connexin 35. Labels gap junctions in retinas of perch, rabbit and macaque in immunohistochemistry.
IMMUNOGEN:	Recombinant Perch Connexin 35
ISOTYPE:	IgG ₁
CLONE NAME:	8F6.2
APPLICATIONS:	Immunohistochemistry Optimal working dilutions must be determined by end user.
SPECIES REACTIVITIES:	Perch, rabbit, and macaque.
FORMAT:	Purified immunoglobulin
PRESENTATION:	Liquid in PBS buffer with 50% glycerol, also contains 0.1% Sodium Azide as a preservative.
STORAGE/HANDLING:	Maintain at -20°C for up to 12 months.

Important Note: During shipment, small volumes of product will occasionally become entrapped in the seal of the product vial. For products with volumes of 200 µL or less, we recommend gently tapping the vial on a hard surface or briefly centrifuging the vial in a tabletop centrifuge to dislodge any liquid in the container's cap.

28820 Single Oak Drive • Temecula, CA 92590
Technical Support: T: 1-800-MILLIPORE (1-800-645-5476) • F: 1-800-437-7502
www.millipore.com



**FOR RESEARCH USE ONLY; NOT FOR USE IN DIAGNOSTIC
PROCEDURES. NOT FOR HUMAN OR ANIMAL CONSUMPTION**

Unless otherwise stated in our catalog or other company documentation accompanying the product(s), our products are intended for research use only and are not to be used for any other purpose, which includes but is not limited to, unauthorized commercial uses, in vitro diagnostic uses, ex vivo or in vivo therapeutic uses or any type of consumption or application to humans or animals.

©2002 - 2012: Millipore Corporation. All rights reserved. No part of these works may be reproduced in any form without permission in writing.