

shelf Life information data sheet

Product Family: Clarigard®

shelf Life information

Shelf Life: 60 Months

Expiration Date: N/A

Storage Conditions: Ambient temperature

Applicable Catalog Number(s): See Attached list

quality Attributes studied

Visual Examination

Water Flow Rate

Gravimetical Extractables

Not all products listed have undergone stability testing to determine the products' shelf lives. The design of the stability schedule is such that only select product(s) within a product family are tested. The design assumes that the products selected are representative of the family as they utilize the same materials of construction, production processes, and packaging components as the products within the family.

The product(s) above does not have a published expiration date or shelf life claim. The shelf life stated above is a guideline based upon available product stability data which is constantly updated as new data is made available. It is not a guarantee or a claim designating the time prior to which a lot of the product will remain within the shelf life specification if stored under the defined conditions, and after which it must not be put into use. It is not an expiration date or a use-by-date. Further, as this is not a product claim, updates to product shelf lives do not require customer notification.

For firms that must establish expiration dates or use-by-dates for this product(s) utilized in a cGMP manufacturing environment, it is recommended that the firms assign a default use-by-date based on their internal procedures/practices.

A handwritten signature in black ink, appearing to read "T. Brackett".

28-Oct-19

Thomas Brackett
Supplier Quality Engineering Manager
Third Party Finished Goods Quality

The life science business of Merck KGaA, Darmstadt, Germany operates as MilliporeSigma in the US and Canada.



D002N0S03	D002N0E03	D002N0V03	D002N0T03	D00271S01	D00271E01
D00271V01	D00271T01	D00201S01	D00201E01	D00201V01	D00201T01
D00251S01	D00251E01	D00251V01	D00251T01	D002F1P01	D00272S01
D00272E01	D00272V01	D00272T01	D00202S01	D00202E01	D00202V01
D00202T01	D00252S01	D00252E01	D00252V01	D00252T01	D002F2P01
D00273S01	D00273E01	D00273V01	D00273T01	D00203S01	D00203E01
D00203V01	D00203T01	D00253S01	D00253E01	D00253V01	D00253T01
D002F3P01	D00274S01	D00274E01	D00274V01	D00274T01	D00204S01
D00204E01	D00204V01	D00204T01	D00254S01	D00254E01	D00254V01
D00254T01	D002F4P01	D003N0S03	D003N0E03	D003N0V03	D003N0T03
D00371S01	D00371E01	D00371V01	D00371T01	D00301S01	D00301E01
D00301V01	D00301T01	D00351S01	D00351E01	D00351V01	D00351T01
D003F1P01	D00372S01	D00372E01	D00372V01	D00372T01	D00302S01
D00302E01	D00302V01	D00302T01	D00352S01	D00352E01	D00352V01
D00352T01	D003F2P01	D00373S01	D00373E01	D00373V01	D00373T01
D00303S01	D00303E01	D00303V01	D00303T01	D00353S01	D00353E01
D00353V01	D00353T01	D003F3P01	D00374S01	D00374E01	D00374V01
D00374T01	D00304S01	D00304E01	D00304V01	D00304T01	D00354S01
D00354E01	D00354V01	D00354T01	D003F4P01	D005N0S03	D005N0E03
D005N0V03	D005N0T03	D00571S01	D00571E01	D00571V01	D00571T01
D00501S01	D00501E01	D00501V01	D00501T01	D00551S01	D00551E01
D00551V01	D00551T01	D005F1P01	D00572S01	D00572E01	D00572V01
D00572T01	D00502S01	D00502E01	D00502V01	D00502T01	D00552S01
D00552E01	D00552V01	D00552T01	D005F2P01	D00573S01	D00573E01
D00573V01	D00573T01	D00503S01	D00503E01	D00503V01	D00503T01
D00553S01	D00553E01	D00553V01	D00553T01	D005F3P01	D00574S01
D00574E01	D00574V01	D00574T01	D00504S01	D00504E01	D00504V01
D00504T01	D00554S01	D00554E01	D00554V01	D00554T01	D005F4P01
D010N0S03	D010N0E03	D010N0V03	D010N0T03	D01071S01	D01071E01
D01071V01	D01071T01	D01001S01	D01001E01	D01001V01	D01001T01
D01051S01	D01051E01	D01051V01	D01051T01	D010F1P01	D01072S01
D01072E01	D01072V01	D01072T01	D01002S01	D01002E01	D01002V01
D01002T01	D01052S01	D01052E01	D01052V01	D01052T01	D010F2P01
D01073S01	D01073E01	D01073V01	D01073T01	D01003S01	D01003E01
D01003V01	D01003T01	D01053S01	D01053E01	D01053V01	D01053T01
D010F3P01	D01074S01	D01074E01	D01074V01	D01074T01	D01004S01
D01004E01	D01004V01	D01004T01	D01054S01	D01054E01	D01054V01
D01054T01	D010F4P01	D030N0S03	D030N0E03	D030N0V03	D030N0T03
D03071S01	D03071E01	D03071V01	D03071T01	D03001S01	D03001E01
D03001V01	D03001T01	D03051S01	D03051E01	D03051V01	D03051T01
D030F1P01	D03072S01	D03072E01	D03072V01	D03072T01	D03002S01
D03002E01	D03002V01	D03002T01	D03052S01	D03052E01	D03052V01
D03052T01	D030F2P01	D03073S01	D03073E01	D03073V01	D03073T01
D03003S01	D03003E01	D03003V01	D03003T01	D03053S01	D03053E01
D03053V01	D03053T01	D030F3P01	D03074S01	D03074E01	D03074V01
D03074T01	D03004S01	D03004E01	D03004V01	D03004T01	D03054S01
D03054E01	D03054V01	D03054T01	D030F4P01		