



3050 Spruce Street
Saint Louis, Missouri 63103 USA
Telephone 800-325-5832 • (314) 771-5765
Fax (314) 286-7828
email: techserv@sial.com
sigma-aldrich.com

Product Information

Deoxyribonucleic acid from herring sperm

Product Number **D 7290**
Storage Temperature -0 °C

Product Description

DNA for hybridization is used to block non-specific binding of probes to membranes.^{1,2} This product is prepared by sonication of purified herring sperm DNA. Sonication shears the large molecular weight DNA to produce fragments in a size range of 587 to 831 base pairs. This range has been shown to be the most effective for hybridizations. The material is monitored during sonication by electrophoresis in order to determine the size range. Once sonication is complete, the material is denatured by boiling.

This is a ready to use concentrated solution (9-12 mg/ml DNA). However, it will reanneal on standing at room temperature so it is recommended to boil the solution for 10 minutes and then cool ON ICE

for at least 5 minutes prior to use. Cooling on ice will reduce the chances for reannealing, as it is more likely to reanneal if cooled at room temperature.

Precautions and Disclaimer

For Laboratory Use Only. Not for drug, household or other uses.

References

1. Sambrook, J., et al., Molecular Cloning: A Laboratory Manual, (Cold Spring Harbor, NY 1989), 9.48-9.49, 13.15.
2. Silhavy, T., et al., Experiments with Gene Fusions, (Cold Spring Harbor, NY: 1984), 195.

MWM/SAG 9/02

Sigma brand products are sold through Sigma-Aldrich, Inc.

Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.