

17137 CPC Selective Supplement

For selective isolation and enumeration of *Vibrio* species from foods.

Composition
(per vial, sufficient for 500 ml medium)

Colistin methane sulfonate	63 mg (680'000 Units)
Polymyxin B sulfate	9 mg (50'000 Units)

Directions

Rehydrate the contents of 1 vial aseptically with 5 ml sterile distilled water. Mix well and aseptically add to 500 ml sterile, molten CPC Agar Base (Prod No. 17134). Mix well and pour into sterile petri plates.

Precautions and Disclaimer

This product is for R&D use only, not for drug, household, or other uses. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

