

## 50279 Listeria UVM II Supplement (Listeria secondary selective enrichment supplement, UVM II)

An antimicrobial supplement recommended for selective cultivation of Listeria species.

### Composition

(per vial, sufficient for 500 ml medium)

Acriflavin hydrchloride 12.5 mg

Nalidixic 10.0 mg

### Directions

Rehydrate the contents of 1 vial aseptically with 2 ml of sterile distilled water. Mix well and aseptically add to 500 ml of sterile, cooled (45-50 °C) UVM Listeria Selective Enrichment Broth, modified I & II. (Cat. No. 90554). Mix well and dispense as desired.

### Precautions and Disclaimer

This product is for R&D use only, not for drug, household, or other uses. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

