

39643 *Burkholderia cepacia* Selective Supplement

Selective supplement for the isolation of *Burkholderia cepacia* in industrial samples

Composition

(per vial sufficient for 500 ml medium)

Polymixin B Sulphate 300,000 IU

Gentamicin 5.00 mg

Vancomycin 1.25 mg

Directions

Reconstitute the vial with 6 ml of sterile distilled water in aseptic conditions and add it to 500 ml of sterilized BCSA Agar base (Cat. No. 39573) cooled to room temperature.

Precautions and Disclaimer

This product is for R&D use only, not for drug, household, or other uses. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

