

P9602 Polymyxin B Selective Supplement

Recommended for selective isolation of various microorganisms.

Composition

(per vial, sufficient for 500 ml or 1000 ml medium)

Polymyxin B Sulphate 50'000 units

Store at 4°C.

Directions

Rehydrate the contents of the vial aseptically with 2 ml of sterile filtered distilled water. Mix well and aseptically add to 500 ml of sterile, molten Cereus Selective Agar (Cat. No 22310), *Bacillus cereus* Agar Base (Cat. No B2176) or to 450 ml of MYP Agar Base (M1928) along with 50 ml Egg Yolk Emulsion (Cat. No 17148). Mix well and pour into sterile petri plates or tubes.

For *Bacillus ChromoSelect* Agar (Cat. No 92325) maximum 1 vial per liter is recommended.

Precautions and Disclaimer

This product is for R&D use only, not for drug, household, or other uses. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

