

18493 MUG Supplement

Fluorogenic substrate for measuring of β -glucuronidase activity, for rapid sensitive identification of *Escherichia coli*.

Composition

(per vial, sufficient for 500 ml medium):

4-Methylumbelliferyl β -D-Glucuronide (MUG) 50.0 mg

Directions

Rehydrate the contents of one vial is 5 ml. of sterile distilled water. Dissolve completely and add to 500 ml of Brilliant Green Bile Lactose Broth (Cat. No 16025) or other mediums. Mix well.

Precautions and Disclaimer

This product is for R&D use only, not for drug, household, or other uses. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

