



QUALITY CONTROL RANGES

ADVANCED BENZONASE CATALOG # EZABNZ-180K STANDARD, LYOPHILIZED CATALOG # E8180-K

STANDARD	Advanced Benzonase® Standard Value
LOT NUMBER: E8180-4K	1000 pg/mL

The exact value of lyophilized standard is shown above.

Reconstitute with **400uL** distilled or deionized water.

QUALITY CONTROLS 1 AND 2, LYOPHILIZED CATALOG #E6180-K

Reconstitute each vial with **250uL** distilled or deionized water.

CONTROL	LOT NUMBER	EXPECTED RANGE	UNITS
QC 1	E6180-104	32 - 67	pg/mL
QC 2	E6180-204	133 - 276	pg/mL

E6180-K: E6180-104.204 (11 December 2025)

[SigmaAldrich.com/support](https://www.sigmaaldrich.com/support) • [SigmaAldrich.com](https://www.sigmaaldrich.com)

© 2025 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved.

The vibrant M logo and MILLIPLEX are trademarks of Merck KGaA, Darmstadt Germany or its affiliates.