

P9352 Perfringens T.S.C. Supplement

Recommended for the selective isolation of *Clostridium perfringens*.

Composition

(per vial, sufficient for 500 ml medium)

Cycloserine 200mg

Store at 4°C.

Directions

Rehydrate the contents of the vial aseptically with 2 ml of sterile filtered distilled water. Mix well and aseptically add it to 475 ml of sterile, molten Perfringens Agar Base (Cat. No 39727) or TSC Agar (Cat. No 93745), along with 25 ml of sterile Egg Yolk Emulsion (Cat. No 17148).

Mix well and pour into sterile petri plates.

Precautions and Disclaimer

This product is for R&D use only, not for drug, household, or other uses. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

