

Product Information

Calcineurin

C1907

Storage Temperature -20 °C

Synonyms

PP2B, Protein Phosphatase 2B, Calcium/Calmodulin-Activated Protein Phosphatase, Calmodulin Binding Protein, Modulator Binding Protein, Phosphoprotein Phosphohydrolase

Product Description

Calcineurin is a cyclosporin-sensitive, calcium- regulated, serine-threonine protein phosphatase with broad substrate specificity. It is the major calmodulin- binding protein found in the brain.¹ This enzyme is a heterodimer composed of catalytic subunit calcineurin A and the regulatory subunit calcineurin B. Calcineurin A has a molecular weight of 58,643 Da and Calcineurin B has molecular weight of 19,200 Da based-on sequence data.^{2,3} Calcineurin A contains four functional domains: a catalytic core with sequence homology to PP-1 and PP-2A, a calmodulin binding site, a calcineurin B binding site and a C-terminal autoinhibitory domain.

Removal of the autoinhibitory domain results in a truncated calcineurin A that is capable of binding the regulatory B subunit and calmodulin, yet no longer requires calcium/calmodulin for full activity.^{4,5}

Calcineurin was first identified as an inhibitor of the calmodulin activation of phosphodiesterase 3':5' cyclic nucleotide (PDE).⁶ This inhibition is caused by the binding of calcineurin to calmodulin which prevents the calmodulin from activating the PDE. Calcineurin has similar effects on adenylate cyclase.⁷ Calcineurin also serves as a key enzyme involved in T-cell activation.^{8,9}

Furthermore, the interaction of cyclosporin A and FK506 with calcineurin is the basis for immuno-suppression by these drugs.^{10,11} Calcineurin is also involved in the hyperphosphorylation of tau protein in Alzheimer's disease¹² and recently has been shown to prevent calpain-mediated proteolysis of tau in differentiated PC12 cells.¹³

Reagent

This is an affinity purified product from bovine brain and supplied as a lyophilized White to light yellow to light brown and faint grey to light grey powder balanced with buffer salts and sucrose as a stabilizer. One unit will cause a 50% inhibition of the activated Phosphodiesterase 3':5' cyclic nucleotide-specific from bovine brain (P9529) activity when assayed with two units of activator (P2277) and 0.1 mM calcium ion in an enzyme coupled system at pH 7.5 at 30 °C.

Precautions and Disclaimer

Please consult the Safety Data Sheet for handling recommendations before working with this material.

Preparation Instructions

The product is soluble in any aqueous buffer.

Storage/Stability

Frost-free freezers are not recommended. We have tested the product by preparing a solution of 50 units per mL in 80 mM Tris, pH 7.5 with 65 mM KCl, 8 mM MgSO₄ and 0.3% albumin. The solution was then stored for seven days at -15 °C then thawed at 4 °C. At the end of the seven days, no percent loss of activity was observed when assayed with phosphodiesterase.

References

1. Klee, C.B., and Krinks, M.H., Biochem., 17, 120- 126 (1978).
2. Ito, A. et al., Biochem. Biophys. Res. Commun., 163, 1492-1497 (1989).
3. Aitken, A. et al., Eur. J. Biochem., 139, 663- 671(1984).
4. Yokoyama, N. et al., Arch. Biochem. Biophys., 300, 615-621 (1993).
5. Hashimoto Y. et al., J. Biol. Chem., 265, 1924-1927 (1990).
6. Wang, J.H. and Desai, R., J. Biol. Chem., 252, 4175-4184 (1977).
7. Antoni F.A. et al., J. Biol. Chem., 270, 28055-28061 (1995).
8. Clipstone, N.A. and Crabtree, G.R., Nature, 357, 695-697 (1992).
9. Batiuk T.D. et al., J Clin. Invest., 100, 1894-1901 (1997).
10. Foor, F. et al., Nature, 360, 682-684 (1992).
11. Liu, J. et al., Biochem. 31, 3896-3901 (1992).
12. Gong, C.X. et al., Brain Res., 741, 95-102 (1996).
13. Xie H.Q. and Johnson G.V, J. Neurosci. Res., 53, 153-164 (1998).

Notice

We provide information and advice to our customers on application technologies and regulatory matters to the best of our knowledge and ability, but without obligation or liability. Existing laws and regulations are to be observed in all cases by our customers. This also applies in respect to any rights of third parties. Our information and advice do not relieve our customers of their own responsibility for checking the suitability of our products for the envisaged purpose.

The information in this document is subject to change without notice and should not be construed as a commitment by the manufacturing or selling entity, or an affiliate. We assume no responsibility for any errors that may appear in this document.

Technical Assistance

Visit the tech service page at SigmaAldrich.com/techservice.

Terms and Conditions of Sale

Warranty, use restrictions, and other conditions of sale may be found at SigmaAldrich.com/terms.

Contact Information

For the location of the office nearest you, go to SigmaAldrich.com/offices.

The life science business of Merck operates as MilliporeSigma in the U.S. and Canada.

Merck and Sigma-Aldrich are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.

© 2024 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved.

C1907pis Rev11/24

