



## QUALITY CONTROL RANGES

### HIGH SENSITIVITY GLUCAGON-LIKE PEPTIDE (GLP-1) ACTIVE CHEMILUMINESCENT ELISA CATALOG # EZGLPHS-35K

#### STANDARD, LYOPHILIZED CATALOG # E8035-HS

Reconstitute standard with **0.5mL** distilled or deionized water.

STANDARD	Active GLP1 Value
LOT NUMBER: 835HS-13K	83 pM

#### QUALITY CONTROLS 1 AND 2, LYOPHILIZED CATALOG # E6035-HS

Reconstitute each vial with **0.5mL** distilled or deionized water.

CONTROL	LOT NUMBER	EXPECTED RANGE	UNITS
QC 1	GLPHS-113	2.3 - 4.8	pM
QC 2	GLPHS-213	11 - 23	pM

**E6035-HS: GLPHS-113.213 (16 December 2025)**

• [SigmaAldrich.com/support](https://www.sigmaldrich.com/support) • [SigmaAldrich.com](https://www.sigmaldrich.com)

© 2025 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved.