


MILLIPORE

Eight-Place Carrier Rack Assembly



**For use with the Millipore
MultiScreen[®] Assay System**

What Is the Eight-Place Carrier Rack Assembly Used For?

The Millipore eight-place carrier rack assembly is a set of 12, eight-place carrier racks in a carrying tray. You use the carrier racks in the Millipore MultiScreen Assay System multiple punch assembly. They hold the sample collection containers in place.

There are three types of eight-place carrier racks available from Millipore; they are for:

- 7 mL vials
- 12 mm X 75 mm test tubes
- 4 mL vials

See the *Millipore MultiScreen Assay System User Guide* for more information on using the MultiScreen system.

Eight-Place Carrier Rack Assembly Catalogue Numbers

Description	Catalogue Number
12, eight-place carrier racks with carrying tray for: 7 mL vials 12 X 75 mm test tubes 4 mL vials	MACR 081 27 MACR 812 75 MACR 081 24

Technical Assistance

Call the office in your country to get additional product information, technical assistance, or order parts.

In the U.S.A., call 1-800-MILLIPORE (800-645-5476). Outside of the U.S.A., see your *Laboratory Catalogue* for the phone number of the Millipore office nearest you.

NOTE: To receive the current edition of our Millipore catalogue, contact Technical Service at the Millipore office closest to you. Or, look up on the Internet at www.millipore.com

Copyright 1997, Millipore Corporation or an affiliated company. Millipore is a registered trademark of Millipore Corporation or an affiliated company. MultiScreen is a trademark of Millipore Corporation or an affiliated company.