



GOAT ANTI-HUMAN APOLIPOPROTEIN A-I POLYCLONAL ANTIBODY

CATALOG NUMBER:	AB740
LOT NUMBER:	
QUANTITY:	1 mL
SPECIFICITY:	Human apolipoprotein A-I. Monospecific by IEP.
IMMUNOGEN:	Purified human apolipoprotein A-I
FORMAT:	Goat serum defibrinated, delipidized and adsorbed by solid phase chromatography as required.
PRESENTATION:	Liquid in 0.05M Tris-HCl, pH 7.5, 0.5M NaCl with 0.1% sodium azide.
STORAGE/HANDLING:	Maintain at -20°C in undiluted aliquots for up to 12 months. Avoid repeated freeze/thaw cycles.
REFERENCE:	Liao, W., et al., <i>Biochemistry</i> (1999) 38 :7532-7544.

Important Note: *During shipment, small volumes of product will occasionally become entrapped in the seal of the product vial. For products with volumes of 200 μ L or less, we recommend gently tapping the vial on a hard surface or briefly centrifuging the vial in a tabletop centrifuge to dislodge any liquid in the container's cap.*

FOR RESEARCH USE ONLY; NOT FOR USE IN DIAGNOSTIC
PROCEDURES. NOT FOR HUMAN OR ANIMAL CONSUMPTION

Unless otherwise stated in our catalog or other company documentation accompanying the product(s), our products are intended for research use only and are not to be used for any other purpose, which includes but is not limited to, unauthorized commercial uses, in vitro diagnostic uses, ex vivo or in vivo therapeutic uses or any type of consumption or application to humans or animals.

©2002 - 2011: Millipore Corporation. All rights reserved. No part of these works may be reproduced in any form without permission in writing.