

61471 ESBL *ChromoSelect* Agar Supplement

ESBL *ChromoSelect* Agar Supplement is used for the detection of Extended Spectrum Beta lactamase producing organisms.

Composition
(per vial sufficient for 500 ml medium)

Ceftazidime	1.50 mg
Cefotaxime	1.50 mg
Ceftriazone	1.00 mg
Aztreonam	1.00 mg
Fluconazole	5.00 mg

Directions:

Rehydrate the contents of 1 vial aseptically with 5 ml sterile distilled water. Mix well and aseptically add to 500ml of sterile, molten, cooled (45-50°C) ESBL *ChromoSelect* Agar (Cat. No. 61471).

Mix well and pour into sterile Petri plates.

Precautions and Disclaimer

This product is for R&D use only, not for drug, household, or other uses. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

