


NovAseptic® Valve Valve Body, Flow through 180° TC

Catalogue Number NA##/32-212 and NU###/32-212

Valve Body, Flow through 180° TC Configuration



Nominal Dimensions in mm (in.)

A	B	C	D*	E	TC	SMS Tube 3008		US Tube ASTM® A269/270		
						NA##	Tube	NU###	Size	Tube
26.8 (1.055)	22.0 (0.866)	42.0 (1.654)	TR40x2	18.5 (0.728)	25.0 (0.984)	NA12	12.0x1.0	NU050	½"	12.7x1.65 (0.50x0.065)
32.0 (1.260)	29.5 (1.161)	52.0 (2.047)	TR55x2	21.0 (0.827)	25.0 (0.984)	NA18	18.0x1.5	NU075	¾"	19.1x1.65 (0.75x0.065)
37.1 (1.461)	39.5 (1.555)	62.0 (2.441)	TR75x2	22.0 (0.866)	50.5 (1.988)	NA25	25.0x1.2	NU100	1"	25.4x1.65 (1.0x0.065)
41.8 (1.646)	52.5 (2.067)	72.0 (2.835)	TR85x2	20.5 (0.807)	50.5 (1.988)	NA38	38.0x1.2	NU150	1½"	38.1x1.65 (1.5x0.065)
59.8 (2.354)	65.5 (2.579)	95.0 (3.740)	TR109x2	32.0 (1.260)	64.0 (2.520)	NA51	51.0x1.2	NU200	2"	50.8x1.65 (2.0x0.065)
75.9 (2.988)	98.0 (3.858)	110.0 (4.331)	TR146x2	34.0 (1.339)	91.0 (3.583)	NA76	76.2x1.65	NU300	3"	76.2x1.65 (3.0x0.065)

* Note: Non-standard thread, special thread to fit NovAseptic Actuators.

Specifications

Net Volume (valve cavity, diaphragm applied)

Valve Body	NA12/ NU050	NA18/ NU075	NA25/ NU100	NA38/ NU150	NA51/ NU200	NA76/ NU300
Volume ml	4	10	30	80	200	470

Net Weight (approximate)

Valve Body	NA12/ NU050	NA18/ NU075	NA25/ NU100	NA38/ NU150	NA51/ NU200	NA76/ NU300
Weight kg (lb)	0.1 (0.22)	0.2 (0.44)	0.5 (1.1)	0.8 (1.8)	1.4 (3.09)	4.3 (9.48)

Material

	Bar Stainless Steel in Compliance with		Tubing Stainless Steel in Compliance with	
Material Code	316L	EN 1.4404	TP316L	EN 1.4404
Technical Requirements	ASME® SA-479	EN 10272	ASTM® A269/270	EN 10217-7/ EN 10216-5

Specifications

Surface Roughness:	Internal surface (manually polished) $Ra \leq 0.5\mu\text{m}$ (20 μin) External surface $Ra \leq 1.6\mu\text{m}$ (63 μin)
Design Temperature, Valve Body:	-80 to 200 °C (-112 to 392 °F)
Design Pressure, Valve Body:	-1 to 10 bar(g) (-14.5 to 145 psi(g))
Note!	The applied diaphragm and actuator may have different design temperature and/or pressure limits. The weakest component in the assembled product determines the maximum design temperature and pressure limits.
Labeling:	Each Valve Body is individually labeled for full traceability and heat No. according to Millipore QA routines.
Packaging:	All open ends on the Valve Body are protected with plastic covers on delivery and packaged in a closed box.
Quality Control:	Millipore Quality Assurance System guarantees the control and traceability at all stages of the manufacturing.
In Compliance With:	The NovAseptic Valve is classified in Europe, under SEP (Sound Engineering Practice) for PED Directive 97/23/EC. NovAseptic Valve is in compliance with ASME® BPE standard.
Options:	For non-standard NovAseptic Valve Body Options, please contact Millipore for further information.

To Place an Order or Receive Technical Assistance

For additional information call your nearest Millipore office or visit us at www.millipore.com.
In the U.S. and Canada, call toll-free 1-800-MILLIPORE (1-800-645-5476). In the U.S., Canada and Puerto Rico, fax orders to 1-800-MILLIFX (1-800-645-5439).

MILLIPORE

Millipore and NovAseptic are registered trademarks of Millipore Corporation or an affiliated company.

ASME is a registered trademark of the American Society of Mechanical Engineers.

ASTM is a registered trademark of the American Society for Testing and Materials.

Specifications subject to change without notice.

Lit. No. SP1131EN00 Rev. 09/2006. Replaces DS-117

© 2006 Millipore Corporation. All rights reserved.