



PURIFIED MOUSE ANTI- HUMAN METAPNEUMOVIRUS (MPV) MONOCLONAL ANTIBODY

CATALOG NUMBER:	MAB8515		
LOT NUMBER:			
QUANTITY:	100 ug	CONCENTRATION:	1 mg/mL
CLONE NAME:	2D7	HOST/ISOTYPE:	IgG ₂
SPECIFICITY:	Reacts against both A and B antigenic groups and subgroups 1 and 2		
IMMUNOGEN:	MPV CAN98-75		
APPLICATIONS:	ELISA IF Optimal dilutions must be determined by the end user		
SPECIES REACTIVITY:	Human Metapneumoviruses		
FORMAT:	0.02 M Phosphate Buffered Saline, 0.25 M Sodium Chloride, 0.1% Sodium Azide, pH = 7.6		
PRESENTATION:	Purified immunoglobulin		
STORAGE/HANDLING:	Store at 2° to 8°C for up to 1 year from date of receipt		
REFERENCES:	Bastien, N. et al. 2003. Sequence analysis of the N, P, M and F genes of Canadian human metapneumovirus strains. <i>Virus Res.</i> May <u>93</u> (1):51-62 Boivin, G. et al. 2002. Virological features and clinical manifestations associated with human metapneumovirus: a new paramyxovirus responsible for acute respiratory-tract infections in all age groups. <i>J. Infect. Dis.</i> <u>186</u> : 1330-1334.		

Important Note: During shipment, small volumes of product will occasionally become entrapped in the seal of the product vial. For products with volumes of 200 µL or less, we recommend gently tapping the vial on a hard surface or briefly centrifuging the vial in a tabletop centrifuge to dislodge any liquid in the container's cap.

FOR RESEARCH USE ONLY; NOT FOR USE IN DIAGNOSTIC
PROCEDURES. NOT FOR HUMAN OR ANIMAL CONSUMPTION

Unless otherwise stated in our catalog or other company documentation accompanying the product(s), our products are intended for research use only and are not to be used for any other purpose, which includes but is not limited to, unauthorized commercial uses, in vitro diagnostic uses, ex vivo or in vivo therapeutic uses or any type of consumption or application to humans or animals.

©2007 - 2009: Millipore Corporation. All rights reserved. No part of these works may be reproduced in any form without permission in writing.

28820 Single Oak Drive • Temecula, CA 92590
Technical Support: T: 1-800-MILLIPORE (1-800-645-5476) • F: 1-800-437-7502
www.millipore.com



upstate | CHEMICON | *Linco*
THE EXPERTISE OF UPSTATE®, CHEMICON® AND LINCO®
IS NOW A PART OF MILLIPORE®
THE EXPERTISE OF UPSTATE®, CHEMICON® AND LINCO®
IS NOW A PART OF MILLIPORE