

3050 Spruce Street, Saint Louis, MO 63103, USA
Website: www.sigmaaldrich.com
Email USA: techserv@sial.com
Outside USA: eurtechserv@sial.com

Product Specification

Product Name:

Product Number: SRE0065

TEST	Specification
Appearance (Turbidity)	Clear
Appearance (Color)	Colorless
Appearance (Form)	Liquid
pH	7.2 - 7.5
at 25°C	
Bioburden	≤ 100 CFU/ml
DNase, Exonuclease Detection	None Detected
RNase Detection	None Detected
NICKase,Endonuclease Detection	None Detected
Protease Detection	None Detected
Osmolality	1640 - 2020 mOs/kg
Endotoxin Level	≤ 1.0 EU/ml

Specification: PRD.0.ZQ5.10000086661

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.