

### Product No. F-2277 Lot 017H0126

# Monoclonal Anti-Mouse TCR α/β FITC Conjugate

Purified Hamster Immunoglobulin Clone H57-597

Monoclonal Anti-Mouse TCR  $\alpha/\beta$  (T Cell Receptor) (hamster IgG) is derived from the H57-597 hybridoma produced by the fusion of P3X63AG8.653 myeloma cells with splenocytes from a hamster immunized with purified murine TCR derived from the DO-11.10 cell line. The product is provided as purified antibody conjugated to fluorescein isothiocyanate isomer I. The conjugate is purified by gel filtration to remove unconjugated FITC. No free FITC is detected. The product is provided at a concentration of 0.25 mg/ml in 0.01 M phosphate buffered saline, pH 7.2, containing 1% BSA, 2 mM EDTA and 0.1% sodium azide (see MSDS)\* as a preservative.

### **Description**

FITC Conjugated Monoclonal Anti-Mouse TCR  $\alpha/\beta$  (clone H57-597) is specific for the 90-110 kD dimer of alpha/beta chains which serve as an antigen specific receptor on T cells in the thymus and periphery. In contrast, the majority of T cells found in epithelial tissue express the gamma/delta chain dimer which is not recognized by H57-597. The TCR dimer is associated with the CD3 complex. In an immobilized form H57-597 will activate T cells to secrete IL-2. Anti-TCR  $\alpha/\beta$  from H57-597 reacts with T cells of all mouse strains tested.

### Performance

When assayed by flow cytometry, approximately 1  $\mu$ g (4  $\mu$ l) of the product will stain 1 x 10<sup>6</sup> mouse spleen cells to maximum fluorescence intensity and percent positive. Prior to adding the product to cells, it is recommended that cell surface Fc receptors be blocked by incubating the cells with 10% - 20% normal goat serum (Sigma Product No. G-9023) for 10 minutes at 4°C. A stabilizing protein such as 1% BSA should be included in the diluent when making dilutions of this product.

**F/P Molar Ratio:** 4.9

### Uses

FITC Conjugated Monoclonal Anti-Mouse TCR  $\alpha/\beta$  may be used for:

- 1. Direct immunofluorescence analysis of mouse T cells by flow cytometry or microscopy.
- Simultaneous multicolor analysis when used in conjunction with PE or Quantum Red<sup>TM</sup> conjugates.

### Storage

Store at 2-8°C. Protect from prolonged exposure to light. If slight turbidity occurs upon prolonged storage, clarify the solution by centrifugation before use.

\* Due to the sodium azide content, a material safety data sheet (MSDS) for this product has been sent to the attention of the safety officer of your institution. Consult the MSDS for information regarding hazardous and safe handling practices.

# Procedure for Direct Immunofluorescent Staining of Splenocytes or Thymocytes

Reagents and Materials Needed but not Supplied

- 1. Diluent: 0.01 M phosphate buffered saline (PBS), pH 7.4, containing 1% BSA and 0.1% NaN<sub>3</sub>.
- 2. FITC conjugated, isotype-matched, non-specific rat or mouse immunoglobulin (Sigma Product No. F-6522).

### Procedure

#### **Notes:**

- 1. In order to obtain best results in different preparations, it is recommended that each individual user determine their optimum dilution of antibody by titration assay.
- 2. Flow cytometric analysis of rodent cells yield better results when the cells are kept cold. Therefore, pre-chill all buffers and diluents, and keep the cells on ice during preparation and staining steps.
- 1. Adjust cell suspension to 1 x 10<sup>7</sup> cells/ml in diluent. Cells should be >90% viable as determined by dye exclusion (e.g., Trypan Blue, Sigma Product No. T-0776). For each sample, add 100 μl or 1 x 10<sup>6</sup> cells per tube. (Note: If cell preparation contains red blood cells, they can be lysed by incubating the cells in approximately 10 mls of 0.017 M Tris, 0.75% NH<sub>4</sub>Cl, pH 7.2 at room temperature for 5 10 minutes followed by centrifugation and washing 2 times in diluent).
- 2. Add 1  $\mu$ g (4  $\mu$ l) of monoclonal antibody to tube(s) containing cells to be stained. Vortex tube gently to mix. Incubate the cells at 4°C for 30 minutes. Proper controls to be included for each sample are:
  - Autofluorescence control: diluent in place of monoclonal antibody, followed by steps 3 - 8.
  - b. Negative staining control: FITC conjugated, isotype-matched, non-specific rat or mouse immunoglobulin at the same concentration as test antibody, followed by steps 3 8.

- After 30 minutes, add 2 ml of cold diluent to all tubes.
- 4. Pellet cells by centrifugation at 500 x G, for 10 minutes.
- 5. Remove supernatant by careful aspiration.
- 6. Resuspend cells in 2 ml of cold diluent.
- 7. Repeat washing procedure (steps 4-6) twice.
- 8. Resuspend cells in 0.5 ml of diluent or 2% paraformaldehyde if cells are to be stored before analyzing.

## **Quality Control**

It is advisable to run the appropriate negative controls. Negative controls establish background fluorescence and non-specific staining of the primary and/or secondary antibodies. The best negative control reagent is a FITC conjugated, isotype-matched, rat or mouse monoclonal antibody or myeloma protein. It should not be reactive with the cells being analyzed and should be used at the same concentration as the fluorophore conjugated specific antibody. The degree of autofluorescence or negative control reagent fluorescence will vary with the type of cells under study and the sensitivity of the instrument used.

### References

Kubo, R., J. Immunol., 142:8, 2736 (1989).

Sigma warrants that its products conform to the information contained in this and other Sigma publications. Purchaser must determine the suitability of the products for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale. Issued 03/97.