



NovaSeptum® Bag Rack for Multi-Sampling Units

Product Description

The NovaSeptum® bag rack for multi-sampling units is designed to easily load and unload NovaSeptum® sampling units. It is fully adjustable and available for all standard holders. A certificate of quality is included in each box.

NovaSeptum® Bag Rack for Multi-Sampling Units



Specifications

Materials of Construction	
Mounting plate	Stainless steel in compliance with ASTM® 316
Sliding plate	Stainless steel in compliance with ASTM® 316
Hanger	Stainless steel in compliance with ASTM® 316
Manifold holder	Polyacetale copolymere
Holder	Stainless steel in compliance with ASTM® 316
Screw	Stainless steel in compliance with A4
Locking nut	Stainless steel in compliance with A4
Wing screw	Stainless steel in compliance with A2/A4

Nominal Dimensions		
Bag rack (length x width x height)	A006	300 mm x 35 mm x 315–515 mm (11.8 in. x 1.3 in. x 12.4–20.2 in.)
	A007	300 mm x 60 mm x 240–440 mm (11.8 in. x 2.3 in. x 9.4–17.3 in.)
	A008	300 mm x 90 mm x 290–490 mm (11.8 in. x 3.5 in. x 11.4–19.3 in.)
	A009	300 mm x 60 mm x 290–490 mm (11.8 in. x 2.3 in. x 11.4–19.3 in.)
	A010	300 mm x 55 mm x 310–510 mm (11.8 in. x 2.1 in. x 12.2–20.0 in.)
	A011	300 mm x 90 mm x 290–490 mm (11.8 in. x 3.5 in. x 11.4–19.3 in.)
	A012	300 mm x 60 mm x 290–490 mm (11.8 in. x 2.3 in. x 11.4–19.3 in.)
Weight	0.6 kg	
Maximum load	5 kg	


Environmental	
Packaging	A package of one bag rack unit with all of its parts is packaged, unassembled, in a cardboard box with a part number.

Ordering Information

NovaSeptum® Bag Rack for Multi-Sampling Units

NovaSeptum® Holder (mm O.D.ø)	Cat. No.
1-port holder, TC 25/NAC 18	A006
1-port holder, TC 25 in-line	A007
5-port holder, TC 50.5 for NAC-A, NAC-D, NAC-I, NAC-S	A008
9-port holder, TC 64 for NAC-A, NAC-D, NAC-I, NAC-S	A009
Ingold® holder	A010
5-port holder, TC 50.5 for NAC-USI, NAC-USM	A011
9-port holder, TC 64 for NAC-USI, NAC-USM	A012

Merck KGaA
Frankfurter Strasse 250
64293 Darmstadt, Germany



To place an order or receive technical assistance

Please visit [MerckMillipore.com/contactPS](https://www.MerckMillipore.com/contactPS)

For additional information, please visit: [MerckMillipore.com](https://www.MerckMillipore.com)